



## **ACTIVE PRODUCT DISPENSING ISOLATOR FOR CHINA**

Extract Technology recently completed the installation and commissioning of another Active Product Dispensing & Mixing Isolator in mainland China, achieving nanogram containment levels.

Xi'an-Janssen Pharmaceutical Ltd were working on a containment project, and had identified the need for an Isolator to protect their operators when handling a new potent compound. The specification and subsequent Purchase Order was based not only on Extract's growing expertise in this area, but also on the success of previous installations for Downflow Booths and Isolators for the customer.

The highly specialised Isolator had many special features, and critical to the success of the project were the ergonomic features that included an inclined front face to the isolator as well as round glove ports and large areas of glazing, all designed with the comfort of the operators in mind and the access they needed.

The Isolator also featured a jacketed vessel with a heater integrated into the base, allowing the operators to dispense the compound straight into the vessel, from where it was transferred under vacuum into a process reactor.

The Xi'an Janssen Project Manager on site commented; "Both the operators and the maintenance personnel are very pleased with the completed installation, and we can confirm that on any future projects that you will be a priority vendor".

***Ends...***

For more information or to contact Extract Technology:

Please call Mike Davies on 01484 317001 or e-mail [mdavies@extract-technology.com](mailto:mdavies@extract-technology.com)

